




*People first*



We believe that it takes more than technology and audiology to create the best hearing instruments. That's why we put the individual needs and wishes of people with hearing loss first in our development of new hearing care solutions.

Oticon ♦ Amigo

Takes you where  
your ears can't go!



With Amigo FM you'll be able to:

- > Defeat noise
- > Zero in on speech
- > Manage your sound environment
- > Be stylish and versatile
- > Communicate with confidence

And all you have to do is point it and listen!

906 55 310 00/05.07 Printed in Denmark



## When noise becomes a problem, Amigo lends a helping hand

Like any other modern individual, you encounter situations that put your ears to the test: a restaurant full of chattering guests, a sidewalk café, where the hum of the crowd and the noise from passing traffic threaten to sabotage every attempt at conversation. Or simply the background noise when watching TV.

When your ears can't capture the sounds you want to hear, Amigo lends a helping hand. This portable FM solution helps you cut through noise so you can get the best out of your hearing wherever you care to go.



### How Amigo works for you

In noisy situations, just two components take you where your ears can't go: The Amigo T10 transmitter, which picks up the person's voice, and the R2 receiver, which plugs into your hearing aid ready to receive a clear speech signal.

### Zero in on speech

- > Amigo's digital sound processor shapes the speech signal for optimal clarity.
- > Soft sounds are boosted and loud sounds controlled, so words sound clearer.

### Defeat noise

- > Amigo uses 4 directional microphones to zero in on the desired speech.
- > Less background noise means that voices stand out more clearly.

### Be stylish and versatile

- > Hi-tech, ergonomic design sits easily in either hand.
- > Supremely versatile – wear it around your neck, carry it or park it on a table.

### Manage your sound environment

- > 3 different operational modes help you cope in the widest variety of situations.
- > Amigo is incredibly easy to use – just turn it on and listen!

### Communicate with confidence

- > High quality components promise excellent long-term durability and reliability.
- > Consistent performance allows you to communicate with confidence.



Open here



# Amigo takes you where your ears can't go!

## Defeat noise

Noise, distance and poor room acoustics can make hearing and understanding a real chore, but the Amigo T10 makes light work of this problem.

Four directional microphones use sophisticated digital sound processing to focus on speech, reject surrounding noise and transmit the signal directly to your hearing aids.

Without the distraction of background sounds you'll experience a dramatic improvement in your ability to follow what family, friends or colleagues are saying.

## Zero in on speech

Team up with the Amigo T10 and be certain that anything said by the people you target will sound closer and so much clearer.

Digital technology gives an extra boost to consonants such as 's' and 'f', that help you to distinguish one word from another. And louder vowel sounds are subdued, to prevent them from smothering the softer speech cues.

Another advantage of digital processing is that the amount of boost given is carefully calculated to match your hearing abilities and listening preferences.



Amigo T10 transmitter and R2 receiver, the portable FM solution



## Be stylish and versatile

You can use the Amigo T10 in conference rooms, classrooms, restaurants and parties, in the street or in your car.

### > Carry it

Held in your hand, the Amigo T10 helps you focus on people in your immediate vicinity.

### > Wear it

Worn around the neck by the person speaking, the Amigo T10 transmits from up to 30m away.

### > Park it

Parked on a table, the Amigo T10 brings people's voices much closer.

### > Connect it

Connect the Amigo T10 to your TV and the whole family will be able to hear at a comfortable volume level.

## Manage your sound environment

The Amigo T10 offers three operational modes to keep you on top of your environment.

When it's noisy, **Focus Mode** is ideal: it captures only the voices appearing directly in its listening zone. This makes it perfect for conversations with two or three friends.

In crowded environments it's time to switch to **SuperFocus**. This mode cuts even more peripheral noise out of the picture, so you can zoom in more easily on a single individual.

In group conversations – at work or at home - **Omni Mode** is perfect as it captures precious details from all directions.

## Communicate with confidence

In the Amigo T10, style and functionality go hand in hand. Ergonomically designed, the shell resembles a high-tech cell phone. Smooth, lightweight and easy to hold, it nestles comfortably in either hand.

The Amigo T10 is built using top quality components and tested to the most rigorous standards, so you can rely on it completely. It will keep you in the loop so that you can communicate with confidence - both now and well into the future.



## Amigo T21 - a professional FM solution

While the Amigo T10 is the finest multi-purpose solution available it's not the only FM solution Oticon offers. If you're taking further education courses or are involved in social or sporting events in venues that could benefit from a truly professional FM system, the Amigo T21 could be the answer.

Easy to use, the Amigo T21 offers all the functionality and reliability you would expect from an Oticon instrument. Your Hearing Care Professional can tell you more about it.